

Refractometer PAL-H



Shop



Quotation Request



Send an inquiry

Mode H™ Technology measure high temperature sample in 2 seconds

Samples with high temperature such as ramen soup, tare sauce for beef bowl, syrup, jam, soy milk etc. had to wait 15 to 20 seconds for measurement. The PAL-H is a high temperature sample, and it has Mode H that can perform stable measurements at a high speed of 2 seconds. "

Specifications

Model	PAL-H
Cat.No.	3870
Range	Brix : 0.0 to 93.0%
Resolution	Brix : 0.5%
Accuracy	±0.5%

Accessories

- MAGIC TM (for PAL MIM) : RE-39446
- STRAP : RE-39410
- PAL-CASE : RE-39409
- JORDAN for PAL(Silicon Cover) : RE-39414
- 10% Sucrose (±0.03%) : RE-110010
- 20% Sucrose (±0.03%) : RE-110020
- 30% Sucrose (±0.03%) : RE-110030
- 40% Sucrose (±0.04%) : RE-110040
- 50% Sucrose (±0.05%) : RE-110050
- 60% Sucrose (±0.05%) : RE-110060
- Battery Case Cover With O-Ring : RE-38163
- AC Battery Charger : RE-34802