



MP300

Automatic Melting Point Apparatus

Hanon MP120 /MP300 Auto Melting Point Apparatus are integrated with photoelectric detection, digital temperature display and other technologies, boasts automatic display of initial and final melts and other functions. In the temperature system, linearly calibrated platinum resistor is used as the detection device, offering efficient and accurate test process and results. The apparatus can be widely used in the chemical industry and medicine research and is essential for the production of medicines, spices, dyes and measurement of other organic crystal substances.

MP300



Maximum Productive
Three independent temperature detection systems ensure the detection of three samples of varied melting points.


Characteristics

- Advanced true-color VGA touch screen is easy for operation and observation.
- 8-level linear temperature rise design is satisfactory for varied temperature rise options.
- Wider temperature control scope and higher temperature accuracy.
- Optimal cooling rate and sample test speed.
- The variation curve of tested sample can be displayed in a timely manner.
- RS232 interface, associated computer software and external computer control are available.
- Detailed test data can be stored.
- Exclusive integral design.

MP120

Automatic Melting Point Apparatus

MP120



- Advanced true-color VGA touch screen is easy for operation and observation.
- 8-level linear temperature rise design is satisfactory for varied temperature rise options.
- Wider temperature control scope and higher temperature accuracy.
- Optimal cooling rate and sample test speed.
- Exclusive integral design.

Note:“●”with the same technical index

Technical data

	MP120	MP300
Melting point range	50℃--400℃	●
"Initial temperature" setting rate	≤6 min(50℃--400℃)	●
	≤7.5 min(400℃--50℃)	●
Temperature resolution	0.1℃	●
Repeatability	0.2℃ at heating rate of 0.2℃/min	●
Temperature accuracy	±0.4℃(<200℃),±0.7℃(<300℃)	●
Temperature gradients	0.1/0.2/0.5/1/2/3/4/5℃/min	0.1-5℃/min 50gradients
Standard capillary size	I.D. 1.0mm, O.D. 1.4mm	●
Capacity per batch	1pc./batch	3pcs./batch
Power supply	220VAC±10% 50Hz	●
Power	150W	300W
Net weight	12kg	16kg
Dimensions	455mmX312mmX200mm	523mmX312mmX200mm